

















T-Port® Enteral Access System - kits and prices

The T-Port® Enteral Access System offers a safe, aesthetic, and well functioning transcutaneous, gastrointestinal access. The system consist of a port made of medical grade titanium and stainless steel and a polyurethane enteral catheter. The design allows easy and rapid placement, and if needed, exchange of catheter. The top of the port, which penetrates the skin, is fitted with a Male Luer-Lock connector allowing easy connection to external devices. The system should be positioned by skilled, and for the purpose, well trained physicians (e.g. interventional radiologist). It is supplied as a sterile device and is intended for single patient use only.

Article	Name (description)	Content			Price (SEK)
1010-1	T-Port® Enteral Access Implant Kit All items required to implant a new system. The components are packed in sterile pouches and delivered with the instructions for use in a single box	T-Port® Base Pack  T-Port® Base	T-Port® Top Pack  T-Port® Top  T-Port® Catheter Cone	T-Port® Tools Pack  Base Holder  Catheter Holder  Top Wrench	20 000
T-Port® Enteral Access Implant Kit T-Port® Enteral Catheter Pack  T-Port® Enteral Catheter					
1020-1	T-port® Enteral Access Replacement Kit All items required to exchange the enteral catheter. The components are packed in sterile pouches and delivered with the instructions for use in a single box	T-Port® Top Pack  T-Port® Top  T-Port® Catheter Cone		T-Port® Tools Pack  Base Holder  Catheter Holder  Top Wrench	10 000
T-Port® Enteral Access Replacement Kit T-Port® Enteral Catheter Pack  T-Port® Enteral Catheter					
1030-1	T-port® Enteral Access Tools Kit All tools required to implant a new system or exchange the enteral catheter. The components are packed in a sterile pouch and delivered with the instructions for use in a single box	T-Port® Enteral Access Tools Kit T-Port® Tools Pack  Base Holder  Catheter Holder  Top Wrench			2 000